

# APC1 Polyclonal Antibody

Cat No: HR1AP2786

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                       |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | APC1 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                              |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                |
| Applications          | WB,IF,ELISA                                                                                                                                                                                                                                                                                                                                           |
| Species Reactivity    | Human,Mouse                                                                                                                                                                                                                                                                                                                                           |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                       |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                       |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                          |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                       |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                       |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                               |
| Observed band         | 215kDa                                                                                                                                                                                                                                                                                                                                                |
| GeneID?Human?         | ANAPC1                                                                                                                                                                                                                                                                                                                                                |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                       |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                       |
| Alternative Names     | ANAPC1; TSG24; Anaphase-promoting complex subunit 1; APC1; Cyclosome subunit 1; Mitotic checkpoint regulator; Testis-specific gene 24 protein                                                                                                                                                                                                         |
| Background            | anaphase promoting complex subunit 1(ANAPC1) Homo sapiens This gene encodes a subunit of the anaphase-promoting complex. This complex is an E3 ubiquitin ligase that regulates progression through the metaphase to anaphase portion of the cell cycle by ubiquitinating proteins which targets them for degradation. [provided by RefSeq, Dec 2011], |