

Spindlin-1 Polyclonal Antibody

Cat No: HR1AP4596

For research use only

Overview

Product Name	Spindlin-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	SPIN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SPIN1; OCR; SPIN; Spindlin-1; Ovarian cancer-related protein
Background	function:May play a role in cell-cycle regulation during the transition from gamete to embryo.,PTM:Phosphorylated during oocyte meiotic maturation.,sequence caution:Translated as Glu.,similarity:Belongs to the SPIN/STSY family.,