



## INTS2 Polyclonal Antibody

Cat No: HR1AP8022

For research use only

### Overview

Product Name	INTS2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	140kDa
GenID?Human?	INTS2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	INTS2; KIAA1287; Integrator complex subunit 2; Int2
Background	integrator complex subunit 2(INTS2) Homo sapiens INTS2 is a subunit of the Integrator complex, which associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690) (Baillat et al., 2005 [PubMed 16239144]).[supplied by OMIM, Mar 2008],