



Cdc20 Polyclonal Antibody

Cat No: HR1AP3035

For research use only

Overview

Product Name	Cdc20 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	CDC20
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDC20; Cell division cycle protein 20 homolog; p55CDC
Background	cell division cycle 20(CDC20) Homo sapiens CDC20 appears to act as a regulatory protein interacting with several other proteins at multiple points in the cell cycle. It is required for two microtubule-dependent processes, nuclear movement prior to anaphase and chromosome separation. [provided by RefSeq, Jul 2008].