

NCAM2 Polyclonal Antibody

Cat No: HR1AP7449

For research use only

Overview

Product Name	NCAM2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	NCAM2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NCAM2; NCAM21; Neural cell adhesion molecule 2; N-CAM-2; NCAM-2
Background	neural cell adhesion molecule 2(NCAM2) Homo sapiens The protein encoded by this gene belongs to the immunoglobulin superfamily. It is a type I membrane protein and may function in selective fasciculation and zone-to-zone projection of the primary olfactory axons. [provided by RefSeq, Jul 2008].