

## Contactin 5 Polyclonal Antibody Cat No: HR1AP7742

For research use only

## Overview

Product Name	Contactin 5 Polyclonal Antibody
Source	Rabbit
Applications	IF,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	••••••
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CNTN5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CNTN5; Contactin-5; Neural recognition molecule NB-2; hNB-2
Background	contactin 5(CNTN5) Homo sapiens The protein encoded by this gene is a member of the immunoglobulin superfamily, and contactin family, which mediate cell surface interactions during nervous system development. This protein is a glycosylphosphatidylinositol (GPI)-anchored neuronal membrane protein that functions as a cell adhesion molecule. It may play a role in the formation of axon connections in the developing nervous system. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2011],