

Contactin 6 Polyclonal Antibody

Cat No: HR1AP6663

For research use only

Overview

Product Name	Contactin 6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CNTN6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CNTN6; Contactin-6; Neural recognition molecule NB-3; hNB-3
Background	contactin 6(CNTN6) Homo sapiens The protein encoded by this gene is a member of the immunoglobulin superfamily. It is a glycosylphosphatidylinositol (GPI)-anchored neuronal membrane protein that functions as a cell adhesion molecule. It may play a role in the formation of axon connections in the developing nervous system. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014],