

## Connexin 46 Polyclonal Antibody

Cat No: HR1AP6638

For research use only

### Overview

Product Name	Connexin 46 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	GJA3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GJA3; Gap junction alpha-3 protein; Connexin-46; Cx46
Background	gap junction protein alpha 3(GJA3) Homo sapiens The protein encoded by this gene is a connexin and is a component of lens fiber gap junctions. Defects in this gene are a cause of zonular pulverulent cataract type 3 (CZP3). [provided by RefSeq, Jan 2010],