

COL9A3 Polyclonal Antibody

Cat No: HR1AP5862

For research use only

Overview

Product Name	COL9A3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	COL9A3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	COL9A3; Collagen alpha-3(IX) chain
Background	collagen type IX alpha 3 chain(COL9A3) Homo sapiens This gene encodes one of the three alpha chains of type IX collagen, the major collagen component of hyaline cartilage. Type IX collagen, a heterotrimeric molecule, is usually found in tissues containing type II collagen, a fibrillar collagen. Mutations in this gene are associated with multiple epiphyseal dysplasia type 3. [provided by RefSeq, Jan 2010],