

COL20A1 Polyclonal Antibody

Cat No: HR1AP8038

For research use only

Overview

Product Name	COL20A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	140kDa
GeneID?Human?	COL20A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	COL20A1; KIAA1510; Collagen alpha-1(XX) chain
Background	COL20A1 (Collagen Type XX Alpha 1) is a Protein Coding gene. Among its related pathways are Collagen biosynthesis and modifying enzymes and ERK Signaling. An important paralog of this gene is MATN1.