



COL23A1 Polyclonal Antibody

Cat No: HR1AP9082

For research use only

Overview

Product Name	COL23A1 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	COL23A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	COL23A1; Collagen alpha-1(XXIII) chain
Background	collagen type XXIII alpha 1 chain(COL23A1) Homo sapiens COL23A1 is a member of the transmembrane collagens, a subfamily of the nonfibrillar collagens that contain a single pass hydrophobic transmembrane domain (Banyard et al., 2003 [PubMed 12644459]).[supplied by OMIM, Mar 2008],