

COL19A1 Polyclonal Antibody

Cat No: HR1AP5870

For research use only

Overview

Product Name	COL19A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	120kDa
GeneID?Human?	COL19A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	COL19A1; Collagen alpha-1(XIX) chain; Collagen alpha-1(Y) chain
Background	collagen type XIX alpha 1 chain(COL19A1) Homo sapiens This gene encodes the alpha chain of type XIX collagen, a member of the FACIT collagen family (fibril-associated collagens with interrupted helices). Although the function of this collagen is not known, other members of this collagen family are found in association with fibril-forming collagens such as type I and II, and serve to maintain the integrity of the extracellular matrix. The transcript produced from this gene has an unusually large 3' UTR which has not been completely sequenced. [provided by RefSeq, Jul 2008],