



CEACAM21 Polyclonal Antibody

Cat No: HR1AP9046

For research use only

Overview

Product Name	CEACAM21 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	CEACAM21
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CEACAM21; Carcinoembryonic antigen-related cell adhesion molecule 21
Background	similarity:Belongs to the immunoglobulin superfamily. CEA family.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,