

Cadherin-24 Polyclonal Antibody

Cat No: HR1AP8280

For research use only

Overview

Product Name	Cadherin-24 Polyclonal Antibody
Source	Rabbit
Applications	IF,WB,ELISA
Species Reactivity	Human, Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CDH24
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDH24; CDH11L; Cadherin-24
Background	function:Cadherins are calcium dependent cell adhesion proteins.,function:Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. Cadherin-24 mediate strong cell-cell adhesion.,similarity:Contains 5 cadherin domains.,subunit:Associates with alpha-, beta- and delta-catenins.,