



Endoglycan Polyclonal Antibody

Cat No: HR1AP9477

For research use only

Overview

Product Name	Endoglycan Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	PODXL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PODXL2; Podocalyxin-like protein 2; Endoglycan
Background	podocalyxin like 2(PODXL2) Homo sapiens This gene is a member of the CD34 family of cell surface transmembrane proteins, which are characterized by an N-terminal extracellular mucin domain, globular and stalk domains, a single pass transmembrane region, and a charged cytoplasmic tail. The encoded protein is a ligand for vascular selectins. [provided by RefSeq, Oct 2012],