

FXR2 Polyclonal Antibody

Cat No: HR1AP9160

For research use only

Overview

Product Name	FXR2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	FXR2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FXR2; FMR1L2; Fragile X mental retardation syndrome-related protein 2
Background	FMR1 autosomal homolog 2(FXR2) Homo sapiens The protein encoded by this gene is a RNA binding protein containing two KH domains and one RCG box, which is similar to FMRP and FXR1. It associates with polyribosomes, predominantly with 60S large ribosomal subunits. This encoded protein may self-associate or interact with FMRP and FXR1. It may have a role in the development of fragile X mental retardation syndrome. [provided by RefSeq, Jul 2008],