

Contactin 1 Polyclonal Antibody

Cat No: HR1AP4963

For research use only

Overview

Product Name	Contactin 1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	113kDa
GeneID?Human?	CNTN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CNTN1; Contactin-1; Glycoprotein gp135; Neural cell surface protein F3
Background	contactin 1(CNTN1) Homo sapiens The protein encoded by this gene is a member of the immunoglobulin superfamily. It is a glycosylphosphatidylinositol (GPI)-anchored neuronal membrane protein that functions as a cell adhesion molecule. It may play a role in the formation of axon connections in the developing nervous system. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011],