

## Nectin 3 Polyclonal Antibody Cat No: HR1AP4023

For research use only

## **Overview**

Product Name	Nectin 3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human, Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	61kDa
GeneID?Human?	PVRL3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PVRL3; PRR3; Poliovirus receptor-related protein 3; CDw113; Nectin-3; CD antigen CD113
Background	nectin cell adhesion molecule 3(NECTIN3) Homo sapiens This gene encodes a member of the nectin family of proteins, which function as adhesion molecules at adherens junctions. This family member interacts with other nectin-like proteins and with afadin, a filamentous actin-binding protein involved in the regulation of directional motility, cell proliferation and survival. This gene plays a role in ocular development involving the ciliary body. Mutations in this gene are believed to result in congenital ocular defects. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011],