



Connexin 25 Polyclonal Antibody

Cat No: HR1AP3138

For research use only

Overview

Product Name	Connexin 25 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	GJB7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GJB7; CX25; Gap junction beta-7 protein; Connexin-25; Cx25
Background	gap junction protein beta 7(GJB7) Homo sapiens Connexins, such as GJB7, are involved in the formation of gap junctions, intercellular conduits that directly connect the cytoplasm of contacting cells. Each gap junction channel is formed by docking of 2 hemichannels, each of which contains 6 connexin subunits (Sohl et al., 2003 [PubMed 12881038]).[supplied by OMIM, Mar 2008],