

## EBFP-Tag Mouse mAb Cat No: HR1AM2285

For research use only

## **Overview**

Product Name	EBFP-Tag Mouse mAb
Source	Mouse
Applications	WB, IHC, IP
Species Reactivity	N/A
Recommended Dilutions	WB 1:5,000 IHC 1:200 IP 1:200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BFP, EBFP, Blue Fluorescent Protein, Enhanced Blue Fluorescent Protein
Background	Enhanced Blue Fluorescent Protein (EBFP) emits a strong blue fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.