

FoxJ3 Polyclonal Antibody

Cat No: HR1AP6379

For research use only

Overview

Product Name	FoxJ3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	FOXJ3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXJ3; KIAA1041; Forkhead box protein J3
Background	similarity:Contains 1 fork-head DNA-binding domain.,