

Paxillin (phospho Tyr31) Polyclonal Antibody

Cat No: HR1AP2504

For research use only

Overview

Product Name	Paxillin (phospho Tyr31) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	68kDa
GeneID?Human?	PXN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PXN; Paxillin
Background	paxillin(PXN) Homo sapiens This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445). [provided by RefSeq, Aug 2011],