

FoxD4/D4L Polyclonal Antibody

Cat No: HR1AP8339

For research use only

Overview

Product Name	FoxD4/D4L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	FOXD4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXD4; FKHL9; FOXD4A; FREAC5; Forkhead box protein D4; Forkhead-related protein FKHL9; Forkhead-related transcription factor 5; FREAC-5; Myeloid factor-alpha; FOXD4L4; FOXD4B; Forkhead box protein D4-
Background	forkhead box D4(FOXD4) Homo sapiens This gene encodes a member of the forkhead/winged helix-box (FOX) family of transcription factors. FOX transcription factors play critical roles in the regulation of multiple processes including metabolism, cell proliferation and gene expression during ontogenesis. Mutations in this gene are associated with a complex phenotype consisting of dilated cardiomyopathy, obsessive-compulsive disorders, and suicidality. [provided by RefSeq, Mar 2012],