

## Vitronectin Polyclonal Antibody Cat No: HR1AP9497

For research use only

## Overview

Product Name	Vitronectin Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	VTN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	VTN; Vitronectin; VN; S-protein; Serum-spreading factor; V75
Background	vitronectin(VTN) Homo sapiens The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond. [provided by RefSeq, Jul 2008],