

## Fibulin-2 Polyclonal Antibody Cat No: HR1AP3454

For research use only

## Overview

Product Name	Fibulin-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	130kDa
GeneID?Human?	FBLN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FBLN2; Fibulin-2; FIBL-2
Background	fibulin 2(FBLN2) Homo sapiens This gene encodes an extracellular matrix protein, which belongs to the fibulin family. This protein binds various extracellular ligands and calcium. It may play a role during organ development, in particular, during the differentiation of heart, skeletal and neuronal structures. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008],