

Integrin β 4 (phospho Tyr1510) Polyclonal Antibody

Cat No: HR1AP7068

For research use only

Overview

Product Name	Integrin β 4 (phospho Tyr1510) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	202kDa
GeneID?Human?	ITGB4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ITGB4; Integrin beta-4; GP150; CD antigen CD104
Background	integrin subunit beta 4(ITGB4) Homo sapiens Integrins are heterodimers comprised of alpha and beta subunits, that are noncovalently associated transmembrane glycoprotein receptors. Different combinations of alpha and beta polypeptides form complexes that vary in their ligand-binding specificities. Integrins mediate cell-matrix or cell-cell adhesion, and transduced signals that regulate gene expression and cell growth. This gene encodes the integrin beta 4 subunit, a receptor for the laminins. This subunit tends to associate with alpha 6 subunit and is likely to play a pivotal role in the biology of invasive carcinoma. Mutations in this gene are associated with epidermolysis bullosa with pyloric atresia. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008],