



Eos Polyclonal Antibody

Cat No: HR1AP8277

For research use only

Overview

Product Name	Eos Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GeneID?Human?	IKZF4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IKZF4; KIAA1782; ZNFN1A4; Zinc finger protein Eos; Ikaros family zinc finger protein 4
Background	IKAROS family zinc finger 4(IKZF4) Homo sapiens Members of the Ikaros (ZNFN1A1; MIM 603023) family of transcription factors, which includes Eos, are expressed in lymphocytes and are implicated in the control of lymphoid development.[supplied by OMIM, Jul 2002],