

ELL2 Polyclonal Antibody

Cat No: HR1AP6387

For research use only

Overview

Product Name	ELL2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	ELL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ELL2; RNA polymerase II elongation factor ELL2
Background	function:Elongation factor that can increase the catalytic rate of RNA polymerase II transcription by suppressing transient pausing by the polymerase at multiple sites along the DNA.,similarity:Belongs to the ELL/occludin family.,subunit:Interacts with EAF1 and EAF2.,