

## COL13A1 Polyclonal Antibody Cat No: HR1AP5867

For research use only

## **Overview**

Decide of Nove	COL 40A4 Palvalanal Antihark
Product Name	COL13A1 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human, Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	COL13A1
Human Swiss- Prot No.	
Cellular localization	
Alternative Names	COL13A1; Collagen alpha-1(XIII) chain; COLXIIIA1
Background	collagen type XIII alpha 1 chain(COL13A1) Homo sapiens This gene encodes the alpha chain of one of the nonfibrillar collagens. The function of this gene product is not known, however, it has been detected at low levels in all connective tissue-producing cells so it may serve a general function in connective tissues. Unlike most of the collagens, which are secreted into the extracellular matrix, collagen XIII contains a transmembrane domain and the protein has been localized to the plasma membrane. The transcripts for this gene undergo complex and extensive splicing involving at least eight exons. Like other collagens, collagen XIII is a trimer; it is not known whether this trimer is composed of one or more than one alpha chain isomer. A number of alternatively spliced transcript variants have been described, but the full length nature of some of them has not been determined. [provided by RefSeq, Jul 2008],