

Cadherin-19 Polyclonal Antibody

Cat No: HR1AP2935

For research use only

Overview

Product Name	Cadherin-19 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	87kDa
GeneID?Human?	CDH19
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDH19; CDH7L2; Cadherin-19
Background	cadherin 19(CDH19) Homo sapiens This gene is one of three related type II cadherin genes situated in a cluster on chromosome 18. The encoded protein is a calcium dependent cell-cell adhesion glycoprotein containing five extracellular cadherin repeats. Loss of cadherins may be associated with cancer formation. Alternative splicing results in multiple transcript variants for this gene. [provided by RefSeq, Aug 2012],