



Nidogen Polyclonal Antibody

Cat No: HR1AP7503

For research use only

Overview

Product Name	Nidogen Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	136kDa
GeneID?Human?	NID1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NID1; NID; Nidogen-1; NID-1; Entactin
Background	nidogen 1(NID1) Homo sapiens This gene encodes a member of the nidogen family of basement membrane glycoproteins. The protein interacts with several other components of basement membranes, and may play a role in cell interactions with the extracellular matrix. [provided by RefSeq, Jul 2008].