

Paxillin Polyclonal Antibody Cat No: HR1AP4267

For research use only

Overview

| Product Name | Paxillin Polyclonal Antibody |
|--------------------------|--|
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 64kDa |
| GeneID?Human? | PXN |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | PXN; Paxillin |
| Background | paxillin(PXN) Homo sapiens This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445). [provided by RefSeq, Aug 2011], |