

COL11A1 Polyclonal Antibody Cat No: HR1AP5863

For research use only

Overview

Product Name	COL11A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human, Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	••••••
Concentration	1 mg/ml
Observed band	181kDa
GeneID?Human?	COL11A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	COL11A1; COLL6; Collagen alpha-1(XI) chain
Background	collagen type XI alpha 1 chain(COL11A1) Homo sapiens This gene encodes one of the two alpha chains of type XI collagen, a minor fibrillar collagen. Type XI collagen is a heterotrimer but the third alpha chain is a post-translationally modified alpha 1 type II chain. Mutations in this gene are associated with type II Stickler syndrome and with Marshall syndrome. A single-nucleotide polymorphism in this gene is also associated with susceptibility to lumbar disc herniation. Multiple transcript variants have been identified for this gene. [provided by RefSeq, Nov 2009],