

Claudin-2 Polyclonal Antibody

Cat No: HR1AP9063

For research use only

Overview

Product Name	Claudin-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	CLDN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CLDN2; PSEC0059; SP82; Claudin-2; SP82
Background	claudin 2(CLDN2) Homo sapiens This gene product belongs to the claudin protein family whose members have been identified as major integral membrane proteins localized exclusively at tight junctions. Claudins are expressed in an organ-specific manner and regulate tissue-specific physiologic properties of tight junctions. This protein is expressed in the intestine. Alternatively spliced transcript variants with different 5' untranslated region have been found for this gene.[provided by RefSeq, Jan 2010].