



COL14A1 Polyclonal Antibody

Cat No: HR1AP8589

For research use only

Overview

Product Name	COL14A1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	COL14A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	COL14A1; UND; Collagen alpha-1(XIV) chain; Undulin
Background	collagen type XIV alpha 1 chain(COL14A1) Homo sapiens This gene encodes the alpha chain of type XIV collagen, a member of the FACIT (fibril-associated collagens with interrupted triple helices) collagen family. Type XIV collagen interacts with the fibril surface and is involved in the regulation of fibrillogenesis. [provided by RefSeq, Jan 2013],