



## DNA pol ? 4 Polyclonal Antibody

Cat No: HR1AP7962

For research use only

### Overview

Product Name	DNA pol ? 4 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GenID?Human?	POLE4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	POLE4; DNA polymerase epsilon subunit 4; DNA polymerase II subunit 4; DNA polymerase epsilon subunit p12
Background	DNA polymerase epsilon 4, accessory subunit(POLE4) Homo sapiens POLE4 is a histone-fold protein that interacts with other histone-fold proteins to bind DNA in a sequence-independent manner. These histone-fold protein dimers combine within larger enzymatic complexes for DNA transcription, replication, and packaging.[supplied by OMIM, Apr 2004],