

CD31 Polyclonal Antibody

Cat No: HR1AP7675

For research use only

Overview

Product Name	CD31 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	130kDa
GeneID?Human?	PECAM1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PECAM1; Platelet endothelial cell adhesion molecule; PECAM-1; EndoCAM; GPIIA ⁺ ; PECA1; CD antigen CD31
Background	platelet and endothelial cell adhesion molecule 1(PECAM1) Homo sapiens The protein encoded by this gene is found on the surface of platelets, monocytes, neutrophils, and some types of T-cells, and makes up a large portion of endothelial cell intercellular junctions. The encoded protein is a member of the immunoglobulin superfamily and is likely involved in leukocyte migration, angiogenesis, and integrin activation. [provided by RefSeq, May 2010],