

COL2A1 Polyclonal Antibody Cat No: HR1AP3133

For research use only

Overview

| Product Name | COL2A1 Polyclonal Antibody |
|--------------------------|---|
| Source | Rabbit |
| Applications | WB,IHC-p,IF,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 140kDa |
| GeneID?Human? | COL2A1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | COL2A1; Collagen alpha-1(II) chain; Alpha-1 type II collagen |
| Background | collagen type II alpha 1 chain(COL2A1) Homo sapiens This gene encodes the alpha-1 chain of type II collagen, a fibrillar collagen found in cartilage and the vitreous humor of the eye. Mutations in this gene are associated with achondrogenesis, chondrodysplasia, early onset familial osteoarthritis, SED congenita, Langer-Saldino achondrogenesis, Kniest dysplasia, Stickler syndrome type I, and spondyloepimetaphyseal dysplasia Strudwick type. In addition, defects in processing chondrocalcin, a calcium binding protein that is the C-propeptide of this collagen molecule, are also associated with chondrodysplasia. There are two transcripts identified for this gene. [provided by RefSeq, Jul 2008], |