

UROC28 Polyclonal Antibody Cat No: HR1AP8127

For research use only

Overview

Product Name	UROC28 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	PBOV1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PBOV1; UROC28; Prostate and breast cancer overexpressed gene 1 protein; Protein UROC28; UC28
Background	prostate and breast cancer overexpressed 1(PBOV1) Homo sapiens This intronless gene encodes a protein of unknown function. Its expression is up-regulated in some types of cancer, including prostate, breast, and bladder cancer. [provided by RefSeq, Aug 2011],