

# UROC28 Polyclonal Antibody

**Cat No: HR1AP8127**

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                    |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | UROC28 Polyclonal Antibody                                                                                                                                                                                                                                         |
| Source                | Rabbit                                                                                                                                                                                                                                                             |
| Applications          | WB,IHC-p,IF,ELISA                                                                                                                                                                                                                                                  |
| Species Reactivity    | Human                                                                                                                                                                                                                                                              |
| Recommended Dilutions |                                                                                                                                                                                                                                                                    |
| Immunogen             |                                                                                                                                                                                                                                                                    |
| Species               | Rabbit                                                                                                                                                                                                                                                             |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                       |
| Isotype               |                                                                                                                                                                                                                                                                    |
| Clonality             |                                                                                                                                                                                                                                                                    |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                            |
| Observed band         | 19kDa                                                                                                                                                                                                                                                              |
| GeneID?Human?         | PBOV1                                                                                                                                                                                                                                                              |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                    |
| Cellular localization |                                                                                                                                                                                                                                                                    |
| Alternative Names     | PBOV1; UROC28; Prostate and breast cancer overexpressed gene 1 protein; Protein UROC28; UC28                                                                                                                                                                       |
| Background            | prostate and breast cancer overexpressed 1(PBOV1) Homo sapiens This intronless gene encodes a protein of unknown function. Its expression is up-regulated in some types of cancer, including prostate, breast, and bladder cancer. [provided by RefSeq, Aug 2011], |