



GFP-Tag Rabbit pAb

Cat No: HR1APEA006

For research use only

Overview

Product Name	GFP-Tag Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	N/A
Recommended Dilutions	WB 1:5,000 IF 1:200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	The green fluorescent protein (GFP) is a protein composed of 238 amino acid residues (26.9kDa) that exhibits bright green fluorescence when exposed to light in the blue to ultraviolet range. In cell and molecular biology, the GFP gene is frequently used as a reporter of expression.