

## EYFP-Tag Mouse pAb Cat No: HR1APEA013

For research use only

## Overview

Product Name	EYFP-Tag Mouse pAb
Source	Mouse
Applications	WB, IP, IF
Species Reactivity	N/A
Recommended Dilutions	WB 1:1,000-3,000 IHC 1:100 IF 1:100
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Yellow florescent protein, Enhanced Yellow florescent protein
Background	Yellow Fluorescent Protein (YFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish Aequorea victoria. Its excitation peak is 514 nm and its emission peak is 527 nm. YFP is a useful tool in cell and molecular biology thanks to its properties useful for fluorescence microscopy.