

?-tubulin Rabbit pAb Cat No: HR1APEA017

For research use only

Overview

Product Name	?-tubulin Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:2,000 IHC 1:200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	52kDakDa
GeneID?Human?	7277
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	Tubulin is one of several members of a small family of globular proteins. The most common members of the tubulin family are ?- tubulin and ?-tubulin, the proteins that make up microtubules. Each has a molecular weight of approximately 55 kiloDaltons. Microtubules are assembled from dimers of ?- and ?-tubulin.