

EFHD1 Rabbit pAb Cat No: HR1APEA021

For research use only

Overview

Product Name	EFHD1 Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Human
Recommended Dilutions	WB 1:1,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG JUUEISU
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	27kDakDa
GenelD?Human?	80303
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MST133, MSTP133, PP3051, Swiprosin 2, SWS2
Background	EFHD1 is an EF-hand domain-containing protein that displays increased expression during neuronal differentiation.