

Glucocorticoid receptor Rabbit pAb Cat No: HR1APEA026

For research use only

Overview

Product Name	Glucocorticoid receptor Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-3,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	80-97kDakDa
GeneID?Human?	2908
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GCCR, GCR, GR, GRL, NR3C1
Background	The glucocorticoid receptor (GR or GCR). The GR is expressed in almost every cell in the body and regulates genes controlling the development, metabolism, and immune response. Because the receptor gene is expressed in several forms, it has many different (pleiotropic) effects in different parts of the body.