



mStrawberry-Tag Rabbit pAb

Cat No: HR1APEA034

For research use only

Overview

Product Name	mStrawberry-Tag Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	N/A
Recommended Dilutions	WB 1:2,000-5,000 IHC 1:200 IF 1:200 IP 1:200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	mStrawberry
Background	mStrawberry is a bright, monomeric red fluorescent protein which was developed in Dr. Roger Tsien's lab by directed mutagenesis of mRFP.