

mOrange-Tag Mouse pAb

Cat No: HR1APEA036

For research use only

Overview

Product Name	mOrange-Tag Mouse pAb
_	
Source	Mouse
Applications	WB
Species Reactivity	N/A
Recommended Dilutions	WB 1:2,000-5,000 IHC 1:200 IF 1:200 IP 1:200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	mOrange
Background	mOrange and mOrange2 are extremely bright orange fluorescent protein monomers which can be used as tags or reporters. Both mOrange fluorescent proteins are mutants derived from mRFP1, a monomeric mutant of DsRed.