

c-Fos Rabbit pAb

Cat No: HR1APEA042

For research use only

Overview

Product Name	c-Fos Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:500-1,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	62kDa
GeneID?Human?	2353
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AP1 antibody, Cellular oncogene c fos antibody, p55 antibody, proto oncogene c fos antibody
Background	c-Fos is a proto-oncogene that is the human homolog of the retroviral oncogene v-fos. c-fos is 380 amino acid protein with a basic leucine zipper region for dimerisation and DNA-binding and a transactivation domain at C-terminus.