

## **? II Tubulin Rabbit pAb** Cat No: HR1APEA048

For research use only

## Overview

Product Name	? II Tubulin Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:2,000-5,000 IHC 1:500-1,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	50kDakDa
GeneID?Human?	7280
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	class II beta tubulin isotype antibody, TUBB 2 antibody, Tubulin beta 2 antibody
Background	Tubulin is the major constituent of microtubules. It exists as a heterodimer consiting of an alpha and a beta subunit. Class II beta tubulin is the major form of tubulin beta in neurons, although it is not neuronal specific. It is found in many other tissues including lung tissue and Schwann cells.