

? III Tubulin Rabbit pAb

Cat No: HR1APEA049

For research use only

Overview

Product Name	? III Tubulin Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:2,000-5,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	50kDaKDa
GeneID?Human?	10381
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	beta 3 tubulin antibody, CDCBM antibody, TUBB3 antibody, Tubulin beta III antibody
Background	Tubulin is the major constituent of microtubules. It exists as a heterodimer consisting of an alpha and a beta subunit. TUBB3 plays a critical role in proper axon guidance and maintenance.